

Order #:	Date:
-----------------	--------------

Name:	Phone #:	Email:
--------------	-----------------	---------------

Bill to:

Account #:	Company name:
-------------------	----------------------

Street:	Street 2:
----------------	------------------

City:	Province:	Postal:
--------------	------------------	----------------

Ship to:

Account #:	Company name:
-------------------	----------------------

Street:	Street 2:
----------------	------------------

City:	Province:	Postal:
--------------	------------------	----------------

Name:	Phone #:	Email:
--------------	-----------------	---------------

Notes:

STANDARD EQUIPMENT (Circle size and choose options below)



- Flex upholstery
- 5 points safety belts with height adjustment
- Seat to back 85° - 95° angle adjustment
- Multi-adjustable foot plates FlipFlap™ with 2 way feet stabilizing belts
- Wheels with glass fiber, front PU swivel, rear inflatable with 2 separate wheel brakes
- Seat - 30° fixed tilt angle
- Backrest
- Adjustable push handles
- Light foldable frame, folding in umbrella shape
- BODYMAP A – seat cushion
- BODYMAP B+ – backrest cushion with side stabilization
- BODYMAP D – headrest
- Reversible electric pump – mobile pump inflating and deflating BODYMAP cushions
- TransFix©








Size 1: \$8,572





Size 2: \$9,100








Size 3: \$9,759








Size 4: \$10,814

ACCESSORIES

CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY
SAFETY	OMO_001 	LED lights - Front and rear LED lights for increased safety.	1	\$119
	OMO_002 	Belt for moving patients AM-P - Robust and ergonomic design as well as high durability facilitates daily care of the patient.	1	\$312
	OMO_003 	Soft cover for tray - Soft cover protects the user from being hit by the tray. Indicated for patients with impaired limbs coordination.	1	\$327
POSITIONING	OMO_107 	Pelvic brace with pelvic belt (adjustable) - Wide belt equipped with "HOLD AND PULL" system. It provides better stabilization in the device.	1	\$234
	OMO_114 	The H Harness acts as a trunk stabilizer. HOLD&PULL© system makes it easier to adjust the vest to the patient body. You can quickly detach the vest if necessary.	1	
		H Harness ('hold & pull' system) (I)		\$186
		H Harness ('hold & pull' system) (A)		\$238
	OMO_115 	Calf belt - Provides additional leg protection. It limits front leg extension and restricts backwards bend	1	\$212
OMO_126 	Trunk belt - Universal solution applicable to many of our products. It secures the user from falling out of the device. Anatomical shape and Velcro straps allow for a comfortable fit.	1	\$256	

	<p>OMO_130</p> 	<p>The 4 points safety vest acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL© system.</p>		
		4 points safety vest (I)	1	\$186
		4 points safety vest (A)		\$256
BODYPAP POSITIONING		<p>*Cotton Covers: In order to minimize the risk of bedsores, we recommend using a cotton cover along with each BodyMap cushion if the permanent duration of use is longer than 2 hours.</p>		
		<p>**VisMemo Covers: For people using BodyMap cushions for more than 8 hours it is recommended to use a thermo-elastic VISMEMO cover.</p>		
		<p>***Aerating-Massaging Pad: It is recommended use of the BodyMap cushion with aerating-massaging pad for people who use BodyMap cushions non-stop longer than 8 hours.</p>		
	OMO_310	U cut Cotton cover forBodymap® A- Sz 2*	1	\$178
		U cut Cotton cover forBodymap® A- Sz 3*		\$178
		U cut Cotton cover forBodymap® A- Sz 4*		\$204
		U cut Cotton cover forBodymap® A- Sz 5*		\$230
	OMO_311	U cut VisMemo cover for Bodymap® A- Sz 2**	1	\$535
		U cut VisMemo cover for Bodymap® A- Sz 3**		\$535
		U cut VisMemo cover for Bodymap® A- Sz 4**		\$609
		U cut VisMemo cover for Bodymap® A- Sz 5**		\$687
	OMO_312	U cut Aerating-massaging pad for Bodymap® A- Sz 2***	1	\$204
		U cut Aerating-massaging pad for Bodymap® A- Sz 3***		\$204
		U cut Aerating-massaging pad for Bodymap® A- Sz 4***		\$230
		U cut Aerating-massaging pad for Bodymap® A- Sz 5***		\$245
	OMO_318	Cotton cover forBodymap® B+ Sz 2*	1	\$223
		Cotton cover forBodymap® B+ Sz 3*		\$223
		Cotton cover forBodymap® B+ Sz 4*		\$249
		Cotton cover forBodymap® B+ Sz 5*		\$282
	OMO_319	VisMemo cover for Bodymap® B+ Sz 2**	1	\$450
		VisMemo cover for Bodymap® B+ Sz 3**		\$450
		VisMemo cover for Bodymap® B+ Sz 4**		\$509
		VisMemo cover for Bodymap® B+ Sz 5**		\$565
	OMO_320	Aerating-massaging pad for Bodymap® B+ Sz 2***	1	\$260
		Aerating-massaging pad for Bodymap® B+ Sz 3***		\$260
		Aerating-massaging pad for Bodymap® B+ Sz 4***		\$290
		Aerating-massaging pad for Bodymap® B+ Sz 5***		\$316
	OMO_322	Cotton cover forBodymap® D – Sz 2	1	\$152
Cotton cover forBodymap® D – Sz 3		\$152		
Cotton cover forBodymap® D – Sz 4		\$230		
Cotton cover forBodymap® D – Sz 5		\$294		

	OMO_323	VisMemo cover for Bodymap® D – Sz 2	1	\$305
		VisMemo cover for Bodymap® D – Sz 3		\$305
		VisMemo cover for Bodymap® D – Sz 4		\$383
		VisMemo cover for Bodymap® D – Sz 5		\$457
COMFORT	OMO_402	 <p>Umbrella - Attachable to the stroller’s frame, perfectly protects from the sunlight or rain.</p>	1	\$375
	OMO_403	 <p>Tray - Easy to assemble/disassemble; can be washed. Element USS_121 is mutually exclusive with USS_405 and USS_420. NOTE: GREY TRAYS ARE AVAILABLE WHILE STOCKS LAST. THEN ONLY BLACK TRAYS WILL BE AVAILABLE.</p>	1	\$553
	OMO_405	 <p>Folding canopy - Protects against sun and rain.</p>	1	\$516
	OMO_408	 <p>Rain cover - Protection from the rain, snow, and wind either for the user and stroller. Element USS_408 and USS_409 available only with USS_405 or USS_420</p>	1	\$398
	OMO_409	 <p>Mosquito net - Protection against insects Element USS_408 and USS_409 available only with USS_405 or USS_420</p>	1	\$93
	OMO_411	 <p>Cup holder - Universal holder mounted to the frame. It enables easy and quick access to bottle or cup.</p>	1	\$163
	OMO_414	 <p>Tray Plexiglass® (2) - Easy to add/remove and can be washed.</p>	1	\$653


	<p>OMO_415</p> 	<p>Rear wheel covers - To cover dirty tires when you want to take the products indoors or transport in trunk.</p>	2	\$26
	<p>OMO_420</p> 	<p>Folding canopy with side covers - Side panels protect against sun, rain and wind.</p>	1	\$632
	<p>OMO_433</p> 	<p>Lifting bar (tip assist) - Makes it easier to tilt the trolley, to go over curbs, steps and other vertical obstacles.</p>	1	\$156
	<p>OMO_439</p> 	<p>Hand brake (4) - The handbrake allows you to stop/slow the buggy with use of hands. It also facilitates maneuvering the trolley. The brake is equipped with a four-stage blockade that also acts as a parking brake.</p>	1	\$669
	<p>OMO_442</p> 	<p>Front wheel covers - Easy to assemble and practical to use. They protect against transferring dirt from tires to indoor spaces.</p>	2	\$19
STORAGE	<p>OMO_502</p> 	<p>Bag DeLux - Large, durable and stylish bag useful for carrying necessary things. Made from high-quality materials.</p>	1	\$702
	<p>OMO_503</p> 	<p>Lady's style bag - This spacious bag designed for mothers to carry all necessities.</p>	1	\$327

	<p>OMO_504</p> 	<p>Man's style bag - Spacious men-style bag for carrying all necessary items.</p>	1	\$327
	<p>OMO_505</p> 	<p>Under seat storage basket - Enables you to carry small objects under stroller. It is attached to the base of the stroller.</p>	1	\$420
	<p>OMO_506</p> 	<p>Carrier bag for the stroller - An excellent solution to protect the folded stroller from damages (e.g. during the transportation). It perfectly fits the stroller shape and dimensions.</p>	1	\$516
	<p>OMO_507</p> 	<p>Oxygen tank holder - Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.</p>	1	Price depends on tank size
	<p>OMO_508</p> 	<p>Saddle bags - Aesthetic and practical side bags mounted at the height of the seat. They allow you to store small items and give the stroller a friendly look.</p>	2	\$375
	<p>OMO_509</p> 	<p>Organizer bag - A handy bag attached to the stroller bar. Specially designed compartments make it easy to keep items in order and the most essential items at hand. The organizer has two openings for a cup or bottle.</p>	1	\$234
BUILD-UP	<p>OMO_601</p> 	<p>Kid-sit Buggypod IO - For convenient and safe transportation of a second child without having to buy a double stroller. The buggypod IO has been designed so that within a few seconds, it can be folded, unfolded, or completely disassembled. 5 point belts seat and folding footrest provide the child with comfort and safety. Using the conversion kit, you are able to convert the kid-sit into an independent three-wheeled stroller. Load capacity up to 20kg.</p>	1	\$1,717

<p>OMO_603</p> 	<p>Extension bar for Buggypod - An adapter enabling to mount the Buggypod kid-sits to the strollers of other manufacturers.</p>	<p>1</p>	<p>\$275</p>
<p>OMO_604</p> 	<p>Mounting set 70 for Buggypod - An adapter enabling to mount the Buggypod kid-sits to the strollers of other manufacturers.</p>	<p>1</p>	<p>\$275</p>
<p>OMO_605</p> 	<p>Rain cover for Buggypod IO - Protection from the rain, snow, and wind either for the user and stroller</p>	<p>1</p>	<p>\$312</p>
<p>OMO_606</p> 	<p>Adaptation set for Buggypod IO - Set of elements which the Buggypod kit-sit can be converted into a fully functional stroller.</p>	<p>1</p>	<p>\$587</p>
<p>OMO_611</p> 	<p>Kid-sit Buggypod PERLE (8) - The award-winning Buggypod PERLE is a perfectly designed seat that fits the BuggyBoard MAXI . When your toddler gets tired, Buggypod PERLE will give him rest and a comfortable ride.</p>	<p>1</p>	<p>\$765</p>
<p>OMO_612</p> 	<p>Buggypod PERLE Adapter - Adapter for mounting Buggypod PERLE to BuggyBoard MAXI .</p>	<p>1</p>	<p>\$230</p>
<p>OMO_613</p> 	<p>Buggyboard MAXI (8) - Useful when the backrest tilt not exceed 60° and seat is facing fwd.</p>	<p>1</p>	<p>\$606</p>

	<p>OMO_622</p>	<p>Buggypod Perle set (6) (8) - The award-winning Buggypod PERLE is a seat that fits the BuggyBoard MAXI . It provides your child with a pleasant ride when tired. The set includes the BuggyBoard MAXI attachment , the Buggypod PERLE seat and the Buggypod PERLE Adapter .</p>	1	\$1,601
SPARE PARTS	<p>OMO_701</p>	Rear inflatable wheel	1	\$186
	<p>OMO_702</p>	Rear inflatable tire	1	\$59
	<p>OMO_703</p>	Rear innertube	1	\$48
	<p>OMO_704</p>	Rear PU wheel (I)	1	\$186
		Rear PU wheel (A)		\$282
	<p>OMO_707</p>	Front inflatable wheel (I)	1	\$186
		Front inflatable wheel (A)		\$282
	<p>OMO_711</p>	Front PU wheel	1	\$282
<p>OMB_736</p>	Adapter 230V	1	\$52	

OMBRELO											
ELEMENT NUMBER	128	129	135	136	405	408	409	412	419	420	421
128											
129											
135											
136											
404											
405						1	1				
408					1					1	
409					1					1	
412											
419											
420						1	1				
421											
ELEMENTS CANNOT BE USED TOGETHER											
ACCESSORY CAN BE MOUNTED UNDER SPECIFIC CONDITIONS											
1. ELEMENTS OMO_408 AND OMO_409 AVAILABLE ONLY WITH CANOPY OMO_405 OR OMO_420											



Symbol	Measurement	Unit	Size 1	Size 2	Size 3	Size 4	Size 5
W	Indicative user height*	[cm]	80 – 100	100 – 130	100 – 150	120 – 180	130 – 180
A/E	Seat/backrest width	[cm]	23	29	34	39	43
L	Seat depth	[cm]	22 – 27	28 – 36	27 – 38	27 – 38	35 – 42
C	Backrest height (with hedrest)	[cm]	43 – 64	46 – 72	46 – 72	53 – 76	54 – 80
F	Footrests length	[cm]	16 – 64	16 – 64	16 – 64	16 – 64	16 – 64
C1	Backrest height (without hedrest)	[cm]	51	58	55	63	60
	Footplate length	[cm]	18	20	21	23	23
	Push handle height	[cm]	100 – 105	100 – 105	100 – 105	100 – 105	100 – 105
X	Backrest tilt angle (relative to the seat)	[°]	90	90	90	90	90
Y	Seat tilt angle	[°]	30	30	30	30	30
	Maximum user weight	[kg]	35	45	55	75	120

Note: Measurement tolerance range is +/- 1 cm / ½ in.

*The value cannot be main factor for selection of product size.



BodyMap A: Seat cushion supports hips in wheelchair, positioning chairs or ordinary chair. What is more, there is also possible to model abduction wedge, raise one hip in order to maintain correct position of spine or off-load sensitive body parts. The shape of seat cushion can be easily changed if

necessary. One side of BM-A cushions is laminated by elastic polyester and the second one is laminated with self-adhesive fabric that allows attachment using additional tapes, straps as well as thigh stabilizers. Self-adhesive function makes it easier to assembly the cushion.



BodyMap B+: B+ backrest is highly recommended for people with back and trunk muscle weakness and for everybody who's in sitting position for a long time, e.g. in the wheelchair. Our positioning cushion reduces excessive stress on skin tissue what limits the risk of bedsores.



BodyMap D Headrest: The orthopedic headrest can be easily shaped to patient's head. It can be mounted to most of the wheelchairs, positioning chairs as well as standing frames. The shape of headrest can be easily changed if necessary. One side of BM-D headrests is laminated by elastic polyester and the second one is laminated with self-adhesive fabric allows to attach additional tapes, straps as well as stabilizers for the patient's head. Self-adhesive function makes it easier to assembly the headrest. If the product is placed in hospital bed, it is necessary to use special bed sheet produced by REH4MAT. This bed sheet has self-adhesive function.