



# RACER EVO BODYMAP special needs stroller

Represented in Canada by: East West Medical Products Ltd.

[www.eastwestmedical.ca](http://www.eastwestmedical.ca)

Email: [info@eastwestmedical.ca](mailto:info@eastwestmedical.ca)

To order / Distributed in Canada by: Assistability

Unit G 7803 – 35 street SE

Calgary AB T2C 1V3

<http://www.assistability.ca/>

Email: [info@assistability.ca](mailto:info@assistability.ca)

Tel: 1-800-403-9019/ 403-629-7697 / Fax: 403.800.9900

<b>Order #:</b>		<b>Date:</b>	
<b>Name:</b>		<b>Phone #:</b>	<b>Email:</b>
<b>Bill to:</b>			
<b>Account #:</b>		<b>Company name:</b>	
<b>Street:</b>		<b>Street 2:</b>	
<b>City:</b>		<b>Province:</b>	<b>Postal:</b>
<b>Ship to:</b>			
<b>Account #:</b>		<b>Company name:</b>	
<b>Street:</b>		<b>Street 2:</b>	
<b>City:</b>		<b>Province:</b>	<b>Postal:</b>
<b>Name:</b>		<b>Phone #:</b>	<b>Email:</b>
<b>Notes:</b>			

## STANDARD EQUIPMENT (Circle size and choose options below)



**Racer Evo**  
(See Racer Evo, Racer, Racer Home order form)



**Racer**  
(See Racer Evo, Racer, Racer Home order form)



**Racer Home**  
(See Racer Evo, Racer, Racer Home order form)



**Racer Evo BodyMap**  
(See p. 2 below)

## Racer Evo BodyMap



- Foldable frame with 4 fixed wheels
- Footrest with feet stabilizing belts
- Seat
- Reclining backrest
- Front rail
- Pump
- Wheels with glass fiber. Front PU swivel with direction blockade, rear inflatable with 2 separate wheel locks
- 1 piece style height adjustable push handle
- 5 points feet stabilizing belts with heel support and foot abduction-adduction adjustments
- 4 points safety vest or pelvic belt
- Lower limbs stabilizer with calf belt
- Leg separator equipped with shin belt
- BODYMAP A+ – seat cushion
- BODYMAP B+ – backrest cushion with side stabilization.
- BODYMAP D – headrest
- Electric pump – mobile pump inflating and deflating BODYMAP cushions
- Lower limb stabilizer





**Size 1:  
\$6,594**










**Size 2:  
\$6,858**

**Size 3:  
\$7,649**









**Size 4:  
\$8,177**

## ACCESSORIES










CATEGORY	SYMBOL	EQUIPMENT	PIECES	QTY
SAFETY	REB_001 	LED lights - Front and rear LED lights for increased safety.	1	\$119
	REB_002 	Belt for moving patients AM-P - Robust and ergonomic design as well as high durability facilitates daily care of the patient.	1	\$312
	REB_003 	Soft cover for tray - Soft cover protects the user from being hit by the tray. Indicated for patients with impaired limbs coordination.	1	\$327
	REB_004 	Low head supporting panels – Sz 1,2 - They provide additional head protection from unwanted contact with stroller's frame. Mounted around the push handle. Element rcr_004 and rcr_005 could not be used together with canopy (el. rcr_404, rcr_405, rcr_420)	1	\$256




		<p>REB_005</p> <p>High head supporting panels, Sz 1,2 - ELEMENT RCR_004 AND RCR_005 COULD NOT BE USED TOGETHER WITH CANOPY (EL. RCR_404, RCR_405 ORAZ RCR_420)</p>	1	\$446
		<p>REB_006</p> <p>Side panels</p>	2	\$539
		<p>REB_107</p> <p>Pelvic brace with pelvic belt (I) - The Pelvic brace with pelvic belt is equipped with "HOLD AND PULL" system. It provides better pelvic stabilization.</p>	1	\$234
		<p>REB_114</p> <p>The H Harness stabilizes the trunk. HOLD&amp;PULL© system makes it easier to adjust the vest to the patient body. Easy to quickly detach the vest.</p>		
		<p>H Harness ('hold &amp; pull' system) (I)</p>	1	\$186
		<p>H Harness ('hold &amp; pull' system) (A)</p>	1	\$238
		<p>REB_126</p> <p>Trunk belt - Universal solution applicable to many of our products. It secures the user from falling out of the device. Anatomical shape and Velcro straps allow for a comfortable fit.</p>	1	\$282
BODYMAP POSITIONING	<p><b>*Cotton Covers:</b> In order to minimize the risk of bedsores, we recommend using a cotton cover along with each BodyMap cushion if the permanent duration of use is longer than 2 hours.</p>			
	<p><b>**VisMemo Covers:</b> For people using BodyMap cushions for more than 8 hours it is recommended to use a thermo-elastic VISMEMO cover.</p>			
	<p><b>***Aerating-Massaging Pad:</b> It is recommended use of the BodyMap cushion with aerating-massaging pad for people who use BodyMap cushions non-stop longer than 8 hours.</p>			
		<p>REB_306*</p> <p>Cotton cover for Bodymap® A+ - Sz 0</p>	1	\$178
		<p>Cotton cover for Bodymap® A+ - Sz 1</p>		\$178
		<p>Cotton cover for Bodymap® A+ - Sz 2</p>		\$204
		<p>Cotton cover for Bodymap® A+ - Sz 3</p>		\$242
		<p>REB_307***</p> <p>VisMemo cover for Bodymap® A+ - Sz 0</p>	1	\$453
		<p>VisMemo cover for Bodymap® A+ - Sz 1</p>		\$502
		<p>VisMemo cover for Bodymap® A+ - Sz 2</p>		\$580
		<p>VisMemo cover for Bodymap® A+ - Sz 3</p>		\$621
		<p>REB_308***</p> <p>Aerating-massaging pad for Bodymap® A+ - Sz 0</p>	1	\$230
		<p>Aerating-massaging pad for Bodymap® A+ - Sz 1</p>		\$230
		<p>Aerating-massaging pad for Bodymap® A+ - Sz 2</p>		\$253
		<p>Aerating-massaging pad for Bodymap® A+ - Sz 3</p>		\$294
		<p>REB_318*</p> <p>Cotton cover for Bodymap® B+ - Sz 0</p>	1	\$175
<p>Cotton cover for Bodymap® B+ - Sz 1</p>		\$189		
<p>Cotton cover for Bodymap® B+ - Sz 2</p>		\$245		
<p>Cotton cover for Bodymap® B+ - Sz 3</p>		\$286		

	<b>REB_319**</b> 	VisMemo cover for Bodymap® B+ - Sz 0	1	\$609
		VisMemo cover for Bodymap® B+ - Sz 1		\$717
		VisMemo cover for Bodymap® B+ - Sz 2		\$873
		VisMemo cover for Bodymap® B+ - Sz 3		\$914
	<b>REB_320***</b> 	Aerating-massaging pad for Bodymap® B+ - Sz 0	1	\$245
		Aerating-massaging pad for Bodymap® B+ - Sz 1		\$282
		Aerating-massaging pad for Bodymap® B+ - Sz 2		\$323
		Aerating-massaging pad for Bodymap® B+ - Sz 3		\$364
	<b>REB_322*</b> 	Cotton cover for Bodymap® D - Sz 0	1	\$119
		Cotton cover for Bodymap® D - Sz 1		\$137
		Cotton cover for Bodymap® D - Sz 2		\$167
		Cotton cover for Bodymap® D - Sz 3		\$208
	<b>REB_323**</b> 	VisMemo cover for Bodymap® D - Sz 0	1	\$279
		VisMemo cover for Bodymap® D - Sz 1		\$301
		VisMemo cover for Bodymap® D - Sz 2		\$346
		VisMemo cover for Bodymap® D - Sz 3		\$386
COMFORT	<b>REB_402</b> 	<b>Umbrella</b> - Attachable to the stroller's frame, perfectly protects from the sunlight or rain.	1	\$375
	<b>REB_403</b> 	<b>Tray - PU</b> - It helps to keep the headrest clean. Easy to assemble/disassemble, it can be washed.	1	\$587
	<b>REB_404</b> 	<b>Folding canopy with leg cover</b> - Protects against both sun and rain.	1	\$843
	<b>REB_405</b> 	<b>Folding canopy</b> - Protects against sun and rain.	1	\$516
	<b>REB_408</b> 	<b>Rain cover</b> - An excellent protection from the rain, snow, and wind for the user and stroller.	1	\$398

	<p>REB_409</p> 	<p><b>Mosquito net</b> - Excellent protection against insects.</p>	1	\$93
	<p>REB_411</p> 	<p><b>Cup holder</b> - Universal holder mounted to the frame. It enables easy and quick access to a drinking container.</p>	1	\$163
Elastic seats and backrests: precisely adapt to the patient's body. Made of memory foam covered with breathable and spatial fabric; Indicated in anti-bedsore prophylaxis.				
	<p>REB_414</p> 	<p><b>Tray Plexiglass®</b> - Easy to add/remove the clear tray can be easily wiped clean.</p>	1	\$743
	<p>REB_415</p> 	For use to cover dirty wheels when transferring stroller indoor or moving in a vehicle.		
		Rear wheel covers (RACER EVO)	2	\$26
		Rear wheel covers (RACER)	4	\$50
	<p>REB_420</p> 	<p><b>Folding canopy with side covers</b> - Side panels protect not only against sun and rain but also the wind.</p>	1	\$632
	<p>REB_442</p> 	<p><b>Front wheel covers (RACER EVO)</b> - Easy to assemble and practical to use. They protect against transferring dirt from tires to indoor spaces.</p>	2	\$19
STORAGE	<p>REB_501</p> 	<p><b>Bag</b> - attaches to stroller bar.</p>	1	\$256
	<p>REB_502</p> 	<p><b>Bag DeLux</b> - Large, durable and stylish bag useful for carrying necessary things. Made from high-quality materials.</p>	1	\$751

	<p>REB_503</p> 	<p><b>Lady's style bag</b> - This spacious bag designed for mothers to carry all necessities.</p>	1	\$327
	<p>REB_504</p> 	<p><b>Man's style bag</b> - Spacious men-style bag for carrying all necessary items.</p>	1	\$327
	<p>REB_505</p> 	<p><b>Under seat storage basket (only RACER EVO)</b> - Enables you to carry small objects under stroller. It is attached to the base of the stroller.</p>	1	\$420
	<p>REB_506</p> 	<p><b>Carrier bag for the stroller</b> - An excellent solution to protect the folded stroller from damages (e.g. during the transportation). It perfectly fits the stroller shape and dimensions.</p>	1	\$557
	<p>REB_507</p> 	<p><b>Oxygen tank holder</b> - Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.</p>	1	Price depends on tank size
	<p>REB_509</p> 	<p><b>Organizer bag</b> - A handy bag attached to the stroller bar. Specially designed compartments make it easy to keep items in order and the most essential items at hand. The organizer has two openings for a cup or bottle.</p>	1	\$234
BUILD-UP	<p>REB_611</p> 	<p><b>Kid-sit Buggypod PERLE (2)</b> - The award-winning Buggypod PERLE is a perfectly designed seat that fits the BuggyBoard MAXI. When your toddler gets tired, Buggypod PERLE will give him rest and a comfortable ride.</p>	1	\$1,115
	<p>REB_612</p> 	<p><b>Buggypod PERLE Adapter</b> - Adapter for mounting <b>Buggypod PERLE</b> to <b>BuggyBoard MAXI</b> .</p>	1	\$230

	<p>REB_613</p> 	<p>Buggyboard MAXI (2) - Useful when the backrest tilt not exceed 70° (size 2) and 60° (size 3)</p>	1	\$606
	<p>REB_622</p> 	<p>Buggypod Perle set (2) (3) - The award-winning Buggypod PERLE is a perfectly designed seat that fits the BuggyBoard MAXI . It provides your child with a pleasant ride when he gets tired.  The set includes the BuggyBoard MAXI attachment , the Buggypod PERLE seat and the Buggypod PERLE Adapter.</p>	1	\$1,784
SPARE PARTS	<p>REB_701</p> 	Rear inflatable wheel (RACER, RACER EVO)	1	\$186
	<p>REB_702</p> 	Rear inflatable tire (RACER, RACER EVO)	1	\$59
	<p>REB_703</p> 	Rear innertube (RACER, RACER EVO)	1	\$48
	<p>REB_704</p> 	Rear PU wheel (I) (RACER, RACER EVO)	1	\$186
		Rear PU wheel (A) (RACER, RACER EVO)	1	\$268
	<p>REB_707</p> 	Front inflatable wheel (RACER EVO) (I)	1	\$186
		Front inflatable wheel (RACER EVO) (A)	1	\$282
	<p>REB_711</p> 	Front PU wheel (RACER EVO)	1	\$282
	<p>REB_719</p> 	Front innertube (RACER, RACER EVO)	1	\$24


	<p>REB_728</p> 	<p>4 points safety vest - It acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&amp;PULL© system.</p>	1	\$245
	<p>REB_736</p> 	Safety rail	1	\$312
	<p>REB_730</p> 	Side tubes protectors	2	\$238
	<p>REB_736</p>	Adapter 230 V	2	\$52

**Attention, some accessories mutually exclude each other. Please refer to the table below.  
Specification of possible collisions between accessories**

<b>RACER / RACER EVO</b>							
ELEMENT NUMBER	004	005	404	405	408	409	420
004							
005							
404					1	1	
405					1	1	
408			1	1			1
409			1	1			1
420					1	1	
<b>ELEMENTS CANNOT BE USED TOGETHER</b>							
<b>ACCESSORY CAN BE MOUNTED UNDER SPECIFIC CONDITIONS</b>							
1. ELEMENTS RCR_408   RCR_409 ARE AVAILABLE ONLY WITH CANOPY (EL. RCR_404, RCR_405 OR RCR_420)							



## SIZE TABLES



Symbol	Measurement	Unit	Size 1	Size 2	Size 3	Size 4
W	Indicative user height*	[cm]	< 120	< 150	< 180	< 180
A	Seat width	[cm]	33	36	39	45
L	Seat depth	[cm]	33	43	47	47
C	Backrest height (with headrest)	[cm]	50 - 62	60 - 72	66 - 84	66 - 84
E	Backrest width	[cm]	30	32	36	41
F	Footrests length	[cm]	30 - 33	31 - 36	37 - 45	37 - 45
C1	Backrest height (without headrest)	[cm]	62	72	84	84
L1	Seat depth to the abduction block	[cm]	15 - 18	23 - 26	28 - 32	28 - 32
	Footer length	[cm]	20	20	22	22
	Push handle height	[cm]	106	106	106	106
Y	Seat tilt angle	[°]	10	10	15	15
X	Backrest tilt angle	[°]	90 - 140	90 - 135	90 - 140	90 - 140
Z	Footplate tilt angle	[°]	100 - 120	100 - 125	100 - 120	100 - 120
	Max. user weight	[kg]	35	50	60	90

**Note:** Measurement tolerance range is +/- 1 cm / ½ in.

\*The value cannot be main factor for selection of product size.



**BodyMap A+ SITTING POSITIONING CUSHION WITH THIGH ABDUCTOR WEDGE AND LATERAL SUPPORTS:** Keeps good posture of patient while sitting. Special abductor wedge reduces contractures of adductors and maintains the legs in the required position. Anatomic shape and lateral supports of BM A+ pillow prevent against the anterior pelvic tilt what provides comfort for bony prominences such as coccyx and sitting bones. This solution reduces the risk of bedsores.



**BodyMap B+** BACKREST CUSHION WITH LATERAL SUPPORTS: highly recommended for people with back and trunk muscles weakness and for everybody who's in sitting position for a long time, e.g. in the wheelchair. Our positioning cushion reduces excessive stress on skin tissue what limits the risk of bedsores.



**BodyMap D HeadRest:** The orthopedic headrest can be easily shaped to patient's head. It can be mounted to most of wheelchairs, positioning chairs as well as standing frames. The shape of headrest can be easily changed if necessary. One side of BM-D headrests is laminated by elastic polyester and the second one is laminated with self-adhesive fabric allows to attach additional tapes, straps as well as lateral supports which stabilize patient's head. Self-adhesive function makes it easier to assembly the headrest.