



Dynamic Arm Life Enablement



The Armon zero gravity Dynamic Arm is designed to help people with limited arm and shoulder function. The products promote free movement of an individual's arms and hands.



The arm support does not take over any function, the user moves themselves. The Armon simply supports the arm, however **unlike a normal armrest, the Armon will dynamically 'follow' the position of the arm through space, providing full support in any position.**

The user only needs some strength in their fingers and some flexibility in their joints: the arm support will eliminate gravity and the user can do the rest.



The Armon dynamic arm supports can assist with:

Training and Rehabilitation - rehabilitation times can be shortened.

Long-term Personal Use - a sophisticated daily-living aide to help the user regain independent function; e.g. eating, drinking, combing hair, putting on glasses, using a phone or a computer.



Dynamic Arm Life Enablement



Examples of conditions which may be helped by Armon Products dynamic arm supports include:

- Muscle disease.
- Neurological disorders like MS, MD, ALS, SMA.
- Brain injury, stroke, CVA, TBI.
- Spinal cord injury.
- Arthritis, rheumatoid arthritis.

We have four different arm supports offering differing levels of functionality:

Mechanical

1. Edero.
2. Pura – steplessly adjustable.

Electrically Driven

3. Elemento.
4. Ayura – tilt mechanism and brakes.

Generally, the less strength the user has, the more advanced the arm support should be. All arm supports can be mounted on a floor stand, chair (including the Bambach Saddle Seat) and wheelchair (manual or electric). The Armon Edero can also be mounted on a table or desk.

Bambach UK

4 Northgate Business Centre, Crown Road, Enfield, London. EN1 1TG

Tel: 020 8532 5100 email: info@bambach.co.uk web: www.bambach.co.uk