

# Pediatric M-FINGER Suggested L-Codes (4 FINGERS and Thumb)



**Partial Hand**  
SOLUTIONS, LLC

<u>Suggested L-Codes</u>	<u>Description</u>
L6020	Partial hand, no fingers remaining
L6680	Test socket, WD/BE
L6687	Frame type socket, WD/BE
L7400	Addition to UE prosthesis, WD/BE Ultralight material (titanium, carbon Fiber or equal)
L7403	Acrylic material, WD/BE
L6695	Addition to upper extremity prosthesis, below elbow/above elbow, custom fabricated from existing mold or prefabricated, socket insert, silicone gel, elastomeric or equal, not for use with locking mechanism <b>***If using a silicone diagnostic socket***</b>
L6696	BE custom fabricated socket insert For congenital or atypical traumatic Amputee, silicone gel, elastomeric or equal For use with or without locking mechanism, Initial use only, use code L6694 or L6695
L7499	M-Finger cable operated mechanical finger, 4 fingers
L7499	Mechanical friction Thumb (M-Thumb)



**PARTIAL HAND SOLUTIONS**  
86 Ladyslipper Lane  
Southington, CT, USA 06489  
(860) 538-5532



**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

