M-FINGERS by Partial Hand SOLUTIONS, LLC



M-Fingers are simple to operate and controlled

by the user's own motion, allowing the patient to perform tasks more effectively, more responsively, and with optimized precision. They conform around objects to provide a secure and confident grip for the user. M-Fingers were designed with high strength, injected molded composite materials to deliver the durability and longevity needed for daily use.







M-FINGER SIZING


Ultra-Small M-Finger







Small Partial M-Finger


















Medium Partial M-Finger



Partial Hand M-FINGER PARTS



M-FINGER PART NUMBER



* (TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: ULTRA-SM-MF-BLK



PARTIAL M-FINGER PART NUMBER

COLOR*

* (TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: SM-PART-BLK





THUMB PART NUMBER

TAN, BRN, BLK

COLOR*

*(TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: THUMB-BLK

2



Partial Hand M-FINGER PARTS SOLUTIONS, LLC

MOUNTING PLATES

	PRODUCT NUMBER	DESCRIPTION
	4F-MP	4 Finger Mounting Plate (in-line)
	3F-MP	3 Finger Mounting Plate (in-line)
	2F-MP	2 Finger Mounting Plate (in-line)
	1F-MP	1 Finger Mounting Plate (in-line)



4 Finger M-Finger System with 4-Finger Mounting Plate, DCM, and BOA mechanism

4 FINGER SYSTEM PART NUMBER

4FS - KIT

COLOR*

TAN, BRN, BLK

SM, MED, LG

* (TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: 4FS-KIT-BLK-SM



4 Finger M-Finger System with 4-Finger Mounting Plate, Thumb, DCM, and BOA mechanism

5 FINGER SYSTEM PART NUMBER

5FS - KIT -

COLOR*
TAN, BRN, BLK

SM, MED, LG

 * (TAN) Tan, (BRN) Brown, or (BLK) Black. Example part number: 5FS-KIT-BLK-SM

HINGES/CONNECTORS

	PRODUCT NUMBER	DESCRIPTION
ii ii	GAFFNEY	Gaffney Joint/Hinge (Colt)
	воа	BOA Fastener for M-Fingers
Q. Q 48	DCM1	1-Finger Direct Connect Mechanism (DCM) Black
	DCM 2	2-Finger Direct Connect Mechanism (DCM) Black
	DCM 3	3-Finger Direct Connect Mechanism (DCM) Black
	DCM 4	4-Finger Direct Connect Mechanism (DCM) Black

GRIP TIPS

PRODUCT NUMBER	DESCRIPTION
GRIPTIP	Silicone protective covering , 5 pack


REPLACEMENT PARTS

REPLACEMEN	II PARIS
PRODUCT NUMBER	DESCRIPTION
PMF-BASE-BLK	Partial M-Finger base (Black)
PMF-BASE-BRN	Partial M-Finger base (Brown)
PMF-BASE-TAN	Partial M-Finger base (Tan)
TB-TIP-BLK	Thumb Tip (Black)
TB-TIP-BRN	Thumb Tip (Brown)
TB-TIP-TAN	Thumb Tip (Tan)
TB-BASE-BLK	Base of thumb (Black)
TB-BASE-BRN	Base of thumb (Brown)
TB-BASE-TAN	Base of thumb (Tan)
SM-MJ-BLK	Small Middle joint (Black)
SM-MJ-BRN	Small Middle Joint (Brown)
SM-MJ-TAN	Small Middle Joint (Tan)
MED-MJ-BLK	Medium Middle joint (Black)
MED-MJ-BRN	Medium Middle Joint (Brown)
MED-MJ-TAN	Medium Middle Joint (Tan)
LG-MJ-BLK	Large Middle joint (Black)
LG-MJ-BRN	Large Middle Joint (Brown)
LG-MJ-TAN	Large Middle Joint (Tan)
FB-BLK	Finger Base Joint (Black)
FB-BRN	Finger Base Joint (Brown)
FB-TAN	Finger Base Joint (Tan)
SM-FT-BLK	Small Finger Tip (Black)
SM-FT-BRN	Small Finger Tip (Brown)
SM-FT-TAN	Small Finger Tip (Tan)
MED-FT-BLK	Medium Finger Tip (Black)
MED-FT-BRN	Medium Finger Tip (Brown)
MED-FT-TAN	Medium Finger Tip (Tan)
LG-FT-BLK	Large Finger Tip (Black)
LG-FT-BRN	Large Finger Tip (Brown)
LG-FT-TAN	Large Finger Tip (Tan)


TITAN SERIES by Partial Hand SOLUTIONS, LLC



The newest innovation from Partial Hand Solutions, the titanium Titan series fingers are designed for the most active heavy duty user, ideal for those with occupational requirements. They are suitable for someone with a full or partial finger amputation and are available in a variety of sizes. The ratchet design allows the patient to manually position the joints in the most functional position, while providing a secure grasp. Various bases allow for the device to be fit as a thumb with full rotation and adjustable friction.



TITAN SIZING







TITAN FINGER PART NUMBER

†(1) Small: length size 3.35", (2) Medium: length size 3.52", or (3) Large: length size 3.62". Example part number: FULL-TITAN3



TITAN PARTIAL FINGER PART NUMBER

PART - TITAN

†(1) Small: length size 2.25", (2) Medium: length size 2.5", or (3) Large: length size 2.62".



TITAN THUMB PART NUMBER

THUMB - TITAN

Example part number: THUMB-TITAN3



M-FINGER ORDER FORM

CUSTOMER #	P.O.#		DATE
PROSTHETIST INFORM.	ATION		
BILLING		SHIPPING (LEAVE BLAN	K IF SAME AS BILLING)
FACILITY/ATTN:		FACILITY/ATTN:	
ADDRESS		ADDRESS	
CITY	STATE/ZIP	CITY	STATE/ZIP
COUNTRY		COUNTRY	
PHONE	FAX	PHONE	FAX
CARRIER* UPS OTH	HER	DATE REQUIRED	TIME
*ORDERS SHIP UPS GROUND C	ON THE FOLLOWING DAY UNLESS SPECIFIED.	EMAIL	
PATIENT INFORMATION	N		
PROSTHETIST NAME		REQUISITIONER	
PATIENT ID / NAME			
NOTES			

For sizing guides, visit college-park.com/fingers





Ultra-Small	ULTRA-SM
Small	SM
Medium	MED
Large	LG
Extra Large	XL

Required Accessories:

- ☐ Gaffney Joint ☐ BOA Fastener
- ☐ Mounting Plate (see below)

MOUNTING PLATES



PARTIAL M-FINGER ALL PARTS INCLUDED






EMAIL COMPLETED ORDER FORM TO: orders@college-park.com **US TOLL FREE** 800.728.7950 **PHONE** 586.294.7950 FAX 586.294.0067 ON THE WEB www.college-park.com (college park

M-FINGER ORDER FORM



For sizing guides, visit college-park.com/fingers





BLK

ADDITIONAL PARTS

DESCRIPTION	PART ID
GRIP TIPS SILICONE PROTECTIVE COVERING, 5 PACK	GRIPTIP
GAFFNEY JOINT/HINGE (COLT)	GAFFNEY
1-FINGER DIRECT CONNECT MECHANISM	DCM1
2-FINGER DIRECT CONNECT MECHANISM	DCM 2
3-FINGER DIRECT CONNECT MECHANISM	DCM 3
4-FINGER DIRECT CONNECT MECHANISM	DCM 4
BOA FASTENER FOR M-FINGER	ВОА

Black




EMAIL COMPLETED ORDER FORM TO: orders@college-park.com **US TOLL FREE** 800.728.7950 **PHONE** 586.294.7950 FAX 586.294.0067 ON THE WEB www.college-park.com



M-FINGER ORDER FORM

PAGE 3

CUSTOMER #	P.O.#		DATE
PROSTHETIST INFORM	ATION		
BILLING		SHIPPING (LEAVE BLAN	K IF SAME AS BILLING)
FACILITY/ATTN:		FACILITY/ATTN:	
ADDRESS		ADDRESS	
CITY	STATE/ZIP	CITY	STATE/ZIP
COUNTRY		COUNTRY	
PHONE	FAX	PHONE	FAX
CARRIER* ☐ UPS ☐ OTI	HER	DATE REQUIRED	TIME
*ORDERS SHIP UPS GROUND ON THE FOLLOWING DAY UNLESS SPECIFIED.		EMAIL	
PATIENT INFORMATIO	N		
PROSTHETIST NAME		REQUISITIONER	
PATIENT ID / NAME			
NOTES			
		For sizing	guides, visit college-park.com/fingers







REPLACEMENT PARTS

PART NUMBER	QUANTITY

PART NUMBER	QUANTITY



EMAIL COMPLETED ORDER FORM TO: orders@college-park.com **US TOLL FREE** 800.728.7950 **PHONE** 586.294.7950 FAX 586.294.0067 ON THE WEB www.college-park.com