

Donning PPE for aerosol generating medical procedures (AGMP) in suspected/confirmed COVID-19 cases


Healthcare professionals must complete these steps prior to performing an aerosol generating medical procedure (AGMP) in patients with suspected or confirmed COVID-19 infection.

Before donning any PPE:

- Remove lanyards, watches, jewelry and other person items. Secure hair back if applicable.
- Perform hand hygiene with alcohol based hand rub. Allow hands to dry before moving to next step.
- Inspect PPE prior to donning. Replace PPE if defects are found.

Recommended PPE for healthcare professionals for AGMPs in a suspected or confirmed COVID-19 case:

- gown
- extended cuff gloves
- N95 respirator (must be fit tested within past 2 years)
- bouffant cap
- face shield



1. Perform hand hygiene



2. Don gown. Secure top and bottom ties



3. Don N95 mask. Ensure straps not crossed



4. Don bouffant



5. Perform N95 mask seal check



5. Don full face shield



5. Don gloves