

Esthetic Multi-Plasma Solution

PLAMERE™



Ultimate various & valuable,

Esthetic Multi-PLASMA Solution

The best assistance for esthetic treatment by plasma technology

Simple and ergonomic design for comfortable grip of stable operation without stress

- All-in-one system with 4 changeable types of plasma treatment probe
- Easy to use by the wheel dial for adjustment of output power strength and range
- Rechargeable battery-powered with micro USB connector



Esthetic Multi-plasma Solution

PLAMERE™


Highlights

Plasma Acne

Acne tip + disposable film

Prevent unnecessary exposure of ozone and UV rays


Treatment by plasma ring with disposable film



Acne tip + acne patch

Plasma ring generates proper amounts of radicals and UV rays under the acne patch

Plasma also helps active ingredients of acne patch to permeate into the skin




Plasma Lifting

Specialized for deep wrinkle

Modulated frequency of 500Hz

Without feedback electrode

HV Induced single-polar probe




Model Name: GPM-1000


Specifications

- Output frequency: 40kHz
- Output power control: analog volume of 1~9 steps
- Mode: continuous, 500Hz modulated pulse
- Various functional tip: lifting ball tip, fractional tip, permeating tip, acne tip with disposable plasma-film
- LED indications: power, shot, low battery, charge status
- Power supply: internal 3.7V 400mAh lithium polymer rechargeable battery with micro USB charger
- Weight: 100g
- Dimension: 165 x 38x 23mm

Plasma Permeating • Fractional




Supports pulse mode



Supports continuous and pulse mode of 500Hz modulated frequency for lifting treatment

Safety power control by analog wheel



Always zero-start for safety after power-on and continuous variable dial of precision control