
PinPointPCR™ No-Extraction Molecular Covid Testing System

Fast. Accurate. Simple.

Turnkey – We provide reagents, equipment, training, validation support

No sourcing needed – ALL items included in per test charge
(swab, plastics, solutions, reagent mixture)

Easy – No extraction. Premixed plates. Sample to result in 1.5 hr.

High-throughput – 752 samples in eight runs on one RT-PCR machine

RT-PCR – The PLATINUM standard for Covid testing

No supply issues – Can produce 100,000 or more assays per order



Saline Swab
(provided)

Mix sample with solution
Heat 5 min



Add heated sample directly to
premixed reaction mixture in 96 well
plate


Load PCR plate into
RT-PCR machine



94 samples per run



1 hour run time
RT-PCR reaction
(not isothermal)



Review data