



ELK Dental Lab


31 O'Neile Crescent NSW 2170
 www.elkdentallab.com
 e-mail: elkdentallab@gmail.com
 Tel: 0480 172 431

Case No.


Work Turnaround Time

Co/Cr frame	10 Days
Teeth Try in	7 Days
Process for Finish	7 Days
Cast Addition	7 Days
Thermoformed	5 Days
Design/Print	5 Days
Weld Point	3 Days
Repair/Reline	2 Day

Please do not include Weekends, Public Holidays or days in transit as working time.

Client:	Tray Base Plate Bite	Date Due:
Patient:		
Cr/Co <input type="checkbox"/> Splint <input type="checkbox"/> Flexible <input type="checkbox"/> Bleaching Tray/s <input type="checkbox"/> Acrylic <input type="checkbox"/> Mouthguard <input type="checkbox"/> Repair/Reline <input type="checkbox"/> Digital/Print <input type="checkbox"/>	Try-In	
Clasp - Design Notation: 	Frame only Frame and Teeth Tooth Setup	
Shade:	Finish Process	

Design Notes:



Proudly manufactured in Australia