



## Simple and effective test method

Used to determine the fastness of a carpet's colour.

Carpet is subjected to a wide variety of testing during everyday wear and tear. To ensure continued high performance colourfastness is critical to the appearance of a carpet. Using the Wira Test Kit, it is possible to determine the fastness of a carpet's colour against standard grey scales and using adjacent fabric.

### Test Method

The test is carried out with a carpet specimen in contact with the specified adjacent fabric under pressure and immersed in water.

The colour change and staining are assessed against grey scales.

The kit is comprised of a perspirometer with perspex plates and weight, oven, 10m multifibre fabric and grey scales.

### Key Features

- ✓ *Perspirometer with plates*
- ✓ *Proven test method*
- ✓ *Oven*
- ✓ *Multifibre Fabric*
- ✓ *Grey Scales*

Standards: ISO 105 E01:1994

Order Code: CCW:001