



Predict cloth shrinkage in making-up with a reproducible measure of 'relaxation' shrinkage

Recommended by leading retailers and manufacturers.

- Provides a rapid routine test for process control
- In extensive use throughout industry and large retailers
- Ergonomic design crafted in fine stainless steel

The Wira Steaming Cylinder subjects four cloth specimens at a time to dry saturated steam at atmospheric pressure. The specimens are given no mechanical restraint and there is no subsequent vacuum cooling. The device provides a reproducible measure of the component of shrinkage occurring in steam pressing attributable to latent strains in the cloth supplied. Any further shrinkage that may occur arises from the characteristics of the pressing process itself.

The Steaming Cylinder is invaluable in maintaining good communications

between clothiers and their suppliers on the specification and performance of woven and knitted fabrics.

Brief Description

The cylinder has a double-walled thermally insulated steaming chamber into which the cloth specimens are placed on a stainless steel supporting frame of low heat capacity and surface area. Steam is fed to the equipment by an integral steam generator. This is connected to a water supply, thus making the operation automatic.

An automatic timer is incorporated in the apparatus to give an indication of the 30-second cycle required for BS4323 and ISO3005. A standard shrinkage rule and coloured pens are supplied to record shrinkage.

Standard Test Method

The specimens are:

- ✓ Preconditioned in the standard atmosphere for testing textiles. Preconditioning ensures that in conditioning itself equilibrium is always approached from the dry side.
- ✓ Marked using the measuring device and then measured.
- ✓ Steamed for 30 seconds, followed by 30 seconds rest outside the cylinder, this cycle being performed three times in all.
- ✓ Preconditioned, conditioned and re-measured.

Conforms to: ISO 3005:1978; BS 4323:1979 [1995]
IWTO-29-2003, Marks & Spencer Standard Test Method P8

Dimensions: 550mm (W), 480mm (D), 340mm (H)

Power Consumption: 3 kW

Consumables:

Replacement sample holder
WIRA Shrinkage rule
Marker pens (set)
Texpens (set)

Order code SCR:SHR
Order code SCR:SRL
Order code SCR:PNS
Order code SCR:TPN

Order Code SCR:001