

ELW Hymns recorded for USB Keyboard Purposes: 170 hymns have been recorded (keyboard) for the purpose of downloading onto a USB stick that can be used on a keyboard. The recorded # is the number that is used for each recording. Two USB sticks will be needed to use all 170 songs, as most keyboards will only allow records 1-99 on each stick. Below are the instructions on how a Clarinova keyboard uses a USB stick with recordings:

#### Music-Keyboard Instructions

Turn keyboard on (right side button)

Push "song select" button on keyboard one time- it should light up "USB."

"S00" should show up on the display. Select the number that you wish to play.

Push Yes to advance to the number you want. Push play to play song. It will stop automatically.

Turn off keyboard following worship.

**YOU CAN NOT RENAME** the songs in the file or on the USB stick. The keyboard only recognizes, for example S00 or S01. Test your recordings prior to using them in worship to make sure they are correct and have downloaded in the correct way.