





- Modular tray design can be scaled to meet your needs.
- Single tray will sanitize 6 masks in several minutes.
- Tray (10" x 17 ½") is small enough to be placed on a standard lab or clinic work surfaces.
- Tray's can be collected and stored on a wheeled rack allowing up to 48 masks to be sanitized at once.
- Heat and chemical free sanitizing. No damage to valve, elastic straps or mask material.
- Ultraviolet C (UVC) light disrupts the virus RNA sequencing killing the virus.
 - Easy to operate. Simple switch and timer function.

For pricing and delivery
Contact ARO Innovation Inc. Toronto Canada 416.824.4716
aroinnovation.com

ALL RIGHTS RESERVED

