

DEBAKEY VASCULAR FORCEPS



Serrated in a longitudinal direction.

Used for grasping soft tissue with minimal damage.

6 1/4" DeBakey Vascular Forceps atraumatic
1.5mm tip
#020-700-16

9 1/2" DeBakey Vascular Forceps atraumatic
1.5mm tip
#17.0624

8" DeBakey Vascular Forceps atraumatic
2mm tip
#17.0820

9 1/2" DeBakey Vascular Forceps atraumatic
2mm tip
#17.0824

6 1/4" DeBakey Vascular Forceps atraumatic
2.7mm tip
#17.0832

12" DeBakey Vascular Forceps atraumatic
2.7mm tip
#17.0835



www.vetsolutionsdirect.com
860-276-7896
866-530-TPLO