




Stool Collection Kit




The Stool Collection Kit is used for collecting and storing fecal samples for various purposes, such as medical diagnostics, research, or laboratory testing. It provides the necessary tools and materials to safely collect, transport, and store fecal specimens for analysis and examination. The collected samples can be used to detect and identify pathogens, parasites, or other substances present in the feces for diagnostic and research purposes.

Storage: Room Temperature

Shelf Life: 3 years

STANDARD COMPONENTS			
NO.	Name	Description	Image
1	Fecal Specimen Catcher	Fecal Specimen Catcher that fits onto any toilet seat, including paper with printed instructions and a plastic wrapping bag	
2	Fecal Collection Tube	Plastic bottle with blue cap and spoon, 30 mL, customized label (Patient information w/ logo)	
3	Cold pack – Single Use	Freezable Gel Pack (3 oz) – single use only	
4	Silver Bubble Mailer	Silver Bubble Mailer Bag (6 x 8 in)	
5	Specimen Bag with Absorbent Pad	Specimen Bag Top Ziplock (6 x 9 in) with Absorbent Pad (3 x 5 in)	
6	Disposable Gloves	Nitrile, latex free, non-sterile, powder-free (3.5 Mil)	
7	Instructions – 2 pages	2 pages – printed Instructions (i.e. Stool Collection Guide, Billing Instructions and etc)	Instructions
8	Packaging Bag (Zip Lock)	Stool Kit Packaging Bag (Zip Lock Type) with customized design (8 x 11 in)	

Packaging Options



Customized box



Customized kit holder for 5 kits



Customized instruction envelope