

EF 6000+ LED LIGHT THERAPY BED

Multi-Wavelength LED Therapy System



The Euro-Fit 6000+ Light Therapy bed is a whole body LED Light Therapy Bed that uses full spectrum of a particular band to treat injuries, reduce inflammation, relax muscle and joints, and increase blood circulation.

Ultimate Full Spectrum 6 Wavelengths:

Advanced Photo biomodulation Therapy Technology:

450nm Blue Light, 520nm Green Light, 633nm Red Light,

(810nm + 850nm + 940nm) Near Infrared Light

HOW DOES LIGHT THERAPY WORK

When delivered within the optimal wavelengths and energy levels, red and near infrared light protect your body's cells from nitric oxide damage, which can otherwise stop the cell's production of ATP when you're stressed or ill. Red light photons allow your cells to continue utilizing oxygen effectively by minimizing the absorption of nitric oxide. Only red light therapy can reach all the way into a cell's mitochondria to stimulate healing and regeneration to help improve your appearance, performance, and overall well-being.

WHY EURO-FIT LED THERAPY LIGHTS ARE SUPERIOR 100+ mW/cm² MOST POWERFUL

The energy output of EURO-FIT 6000+ LED Therapy Bed lights are maximized, ensuring you get a clinical-grade red light therapy treatment.

Our LED lights emit more irradiance (the output of an LED therapy light) than most of the other LED Therapy Lights currently on the market.

In addition to superior output strength, EURO-FIT LEDs offer a highly advanced Red & Near Infra-Red wave spectrum comprised of four independent and highly studied bands of light.

LED PROXIMITY:

Distance from the LED affects light absorption, but do not be fooled by erroneous statements of needing to be directly on the LEDs. That is absolutely necessary for weaker beds that do not use 220v, high irradiance and a high-density of LEDs.

LIFETIME LEDs:

EURO-FIT 6000+ utilizes 100,000 hour, German LEDs that are under warranty for life. That's right. You're covered for life if any LEDs go out.

Note: 100K hour LEDs are the max lifespan in the industry. Claims of higher lifespans are unproven.

EURO-FIT 6000+ LED THERAPY BED HIGHLIGHTS

Wireless Control Panel:

The EURO-FIT 6000+ LED Light Therapy bed offers adjustability of irradiance and frequency for each of the six wavelengths independently. All done with a wireless hand-held tablet.

360 Degree Whole Body Design:

The XL surface allows even a 7ft. Person to lie in it without obstacles and the professional design and the right distance to make sure the Red/NIR energy to be absorbed efficiently and evenly.

Surround Sound Plus:

The EURO-FIT 6000+ LED Light Therapy Bed comes with pre-installed speakers which has subwoofer function and 2 extra music channels. Couple this with MP3 integration and users will have an audio experience they will love.

BENEFITS TO USING RED LIGHT THERAPY LIGHT

曲

STIMULATE HAIR GROWTH

 Stimulates stem cells in the hair follicle to encourage active growth

IMPROVED SKIN HEALTH

- · Generates production of collagen
- Increased tissue repair
- Reduction in appearance of fine lines, wrinkles, and scars

PROMOTE CELLUAR HEALTH

- Stimulates collagen production
 Strengthens hair and improves skin elasticity
- Promotes connective tissue health



IMPROVE CIRCULATION

- Dilates blood vessels for increased flow
- Protects red blood cells and platelets

REDUCE RECOVERY TIME

- · Accelerates muscle repair
- Stimulates mitochondria and stem cells for faster recovery

REDUCE INFLAMMATION


- Soothes sore muscles, joint pain, and arthritis
- Eases symptoms associated with autoimmune diseases, spinal cord and traumatic brain injuries

INCREASE FERTILITY

Increases testosterone production
May stimulate the brain's pineal gland

REDUCE PAIN

- Eases joint stiffness and soreness
- Diminishes inflammation
- Reduces muscle spasms and increases blood flow



HARNESS THE POWER OF RED | NIR

PRODUCT PARAMETERS

MODEL	EURO-FIT EF 6000+ LED Light Therapy Bed
LIGHT TYPE	41600 LEDs Made in Germany
LED RATING	100,000 HRS
WAVELENGTHS UTILIZED	450nm, 520nm, 630nm, 810nm, 850nm, 940nm Ability to independently control each wavelength when needed
IRRADIANCE	120mW/cm ² Adjustable control 1-120W/cm ²
RECOMMENDED TREATMENT TIME	10-20 min
TOTAL DOSE IN 10 MIN	60J/cm ²
OPERATION MODE	True Continuous Wave Variable Pulse 1-5000Hz in 1Hz Increments Ability to Change Pulse
WIRELESS TABLET CONTROL	Ability to Manage Multiple Systems Ability to Set and Store Protocols Ability to Control from Front Desk
PRODUCT SPECIFICATIONS	Dimensions: 2198mm X 1157mm X 1079mm (closed) Weight: 300Kg Weight Capacity: 300 Kg
ELECTRICAL REQUIREMENTS	220-240 VAC 50/60Hz 30A Single Phase
FEATURES	360 Degree Treatment Reflective Panels Homogeneous Light Distribution Air Cooling System Bottom Wheels for Mobility Built-in Bluetooh Speakers
WARRANTY	2 Years



www.geteuro-fit.com

832-894-0580