



Phone: 203-433-4329  
Fax: 203-208-0664  
Toll Free: 866-702-0260  
[www.phoenixsongsbio.com](http://www.phoenixsongsbio.com)

## **Cryopreserved Hepatocyte Recovery Medium**

### **Complete Medium for the Recovery of Mature Cryopreserved Human Hepatocytes**

Cryopreserved Hepatocyte Recovery Medium (CRM) is formulated for the optimal recovery of adult and neonatal cryopreserved human hepatocytes. CRM is shipped frozen on dry ice and upon thawing is complete and ready to use for recovery of cryopreserved hepatocytes. . Thaw 1-2 vials of cryopreserved human hepatocytes per 50ml of CRM.

### **Materials Needed When Thawing Cryopreserved Human Hepatocytes**

1. Cryopreserved Hepatocyte Recovery Medium (CRM050)
  - a. Can thaw 1-2 vials of hepatocytes into 50ml of CRM050
  - b. CRM050 comes frozen and is ready for use upon thawing
2. Cryopreserved Hepatocyte Plating Medium (HPM250)
  - a. Make according to the media insert accompanying the medium
3. Hepatocyte Culture Medium (Cat.# HCM500)
  - a. Make according to the media insert accompanying the medium
4. Follow the Recovery and Culture of Cryopreserved Hepatocytes Methods for full details.

**Use Restrictions:** The Products and components of the Products are to be used for RESEARCH PURPOSES only and are not for use in humans or animals. No rights expressed or implied are conveyed by the sale of the Products to make, have made, offer to sell or sell the Products alone or in combination with any other products or to use for commercial purposes without prior written consent.