

T-Edge Autoclaves

10" and 11" Chamber Sizes

Gravity autoclave upgradable
to pre/post vacuum

Maximum efficiency, Minimum costs.

- **35 Minute wrapped cycle time including drying!***
- Optimal chamber capacity (more than 80%), with extra-long trays
- Convenient and easy loading with door opening to 110°
- User-friendly touchpad and control panel for displaying programs, cycles times, water levels and more
- Advanced documentation and traceability options and barcode
- Remote Wi-Fi operating via smartphone or tablet

**Cycle time based on T-Edge 10*



ONLY
UPGRADABLE
AUTOCLAVE
IN THE
MARKET!

Designed to meet the most current
sterilization standards ANSI/AAMI ST55


Tuttnauer
Innovation · Legacy · Partnership

Advanced barcode traceability



Receive barcodes with sterilization data (including expiration date).



Apply barcode label to each pouch.



Scan barcode data for patient file during treatment, for safety and regulatory purposes.

Sterilization Programs

Icon	Cycle Program	Load Type	Cycle Time (min.)*	Drying Time (min.)
	Unwrapped 134	Unwrapped instruments	T-10: 25 min T-11: N/A	T-10: 1 min T-11: N/A
	Wrapped 134	Wrapped instruments	T-10: 25 min T-11: N/A	T-10: 10 min T-11: N/A
	Unwrapped 121	Unwrapped delicate instruments	T-10: 39 min T-11: N/A	T-10: 1 min T-11: N/A
	Wrapped 121	Wrapped delicate instruments	T-10: 50 min T-11: N/A	T-10: 10 min T-11: N/A

* Cycle times are based on: 1 kg standard loads (not including drying), Warm chamber (and subject to load size)

Technical Specifications

Autoclave Model	Chamber Size (cm)	Volume	Overall Dimensions D x W x H (cm)	Number of trays	Connections	Item Supplies	Power Supply	Power Consumption	Weight
T-Edge 10	18" depth x 10" diameter	6 gal/ 23 liters	18.90"x19.23" x22.44"	5	1 ethernet port, 4 USA ports, WiFi, integrated automatic water filling and drain connection	Drain tube, tray holder	230v, Frq 60Hz	2000W	110 lbs.
T-Edge 11	18" depth x 11" diameter	7 gal/ 27 liters	18.90"x19.69" x22.44"	5	1 ethernet port, 4 USA ports, WiFi, integrated automatic water filling and drain connection	Drain tube, tray holder	120V, Frq 60Hz	1500W	115 lbs.



Learn from our experts.
Visit our blog: tuttnauer.com/blog



Tuttnauer USA Co. Ltd.
25 Power Drive, Hauppauge, NY 11788
Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720
E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com

