



**THE SS-24 UNIT UTILIZES A TUNED RESONATE CHAMBER COMPRISED OF AN ARRAY OF DI-POLAR MATERIALS IN A FLOW FORM. AS THE WORT MOVES THROUGH THE CHAMBER IN A DOUBLE HELIX MOTION, A SERIES OF TEMPERATURE AND PRESSURE FLUCTUATIONS OCCUR WITHIN THE WORT WHICH ORGANIZE PARTICLES AND CREATES A BIO-MAGNETIC FIELD IN AND AROUND THE WORT FORMING AN OPTIMAL ENVIRONMENT FOR HEALTHY YEAST CELLS. THE WORT HAS MORE AVAILABLE OXYGEN AND ENHANCES NUTRIENT ABSORPTION FOR GENERATIONAL YEAST VIABILITY.**

