

FEATURED PRODUCT



Edition: May 2026

As social safety concerns rise, our portable self-protection devices make beverage screening fast and easy.

Ideal for on-the-go use, our wristbands require no additional equipment. Simply apply one or three drops of your drink, and obtain clear, accurate results in just 3 minutes through intuitive visual cues. Featuring precise cut-off values and high detection accuracy, they reliably identify drugs including GHB, KET, FLZ, SCOP and COC.

More than just a testing tool, it puts personal safety in your hands, supporting quick risk identification and enhancing safety awareness in any social environment.

PRODUCT ADVANTAGES

- ✓ **Rapid Detection** : Delivers results within 1–3 minutes via a simple color change reaction.
- ✓ **User-Friendly** : 3-step operation requires no technical skills.
- ✓ **Portable** : Compact, lightweight, and easy to carry in daily life or travel.
- ✓ **Multi-Substance Detection** : Screening for GHB, KET, FLZ, COC and SCOP

PRODUCT INFORMATION

DrinkSafe Check - GHB&KET

Cat. No.: D26-MDT (B)

Specimen: Drink

Test Size: 10 T

Cut-off: GHB: 12mg/mL,

KET: 4mg/mL

DrinkSafe Check - FLZ&GHB&KET

Cat. No.: D36-MDT (B)

Specimen: Drink

Test Size: 10 T

Cut-off: FLZ: 250µg/mL,

GHB:12mg/mL, KET:4mg/mL

DrinkSafe Check - SCOP&FLZ&GHB&KET/COC

Cat. No.: D56-MDT (B)

Specimen: Drink

Test Size: 10 T

Cut-off: SCOP:2000ng/mL,

FLZ: 250µg/mL, GHB:12mg/mL, KET/COC:4mg/mL



Featured Product Recommendation

DRINKSAFE CHECK WRISTBANDS