

Caremed Alrick.

Mobilise® CH5-L bariatric chair

An innovative early mobilisation chair*, made in the UK



An innovative bariatric, early mobilisation chair that lays 100% flat and has an innovative legrest that encourages sit to standing, this chair is manufactured in the UK, using UK materials which greatly reduces leadtimes.

Manual Handling

100% Layflat® reduces the need for hoists



Easy to move

Allows patients to be moved outside or to other wards if necessary



Seating options

Depending on the needs of the patient, different seating options are available



Retractable legrest

Fully retracts beneath the seat so easier sit to standing



Bariatric version

Seat width to 660mm and safe working load of 250kg (39 stone)



Made in the UK

The Mobilise[®]chair is made in the UK



Features

- Tilt in space chair
- ✓ Lays 100% flat
- ✓ Sit to standing option
- Retractable legrest fully retracts beneath the seat so easier sit to standing
- Seat tilt with Trendelenburg
- Easily remove the seat and back cushions for cleaning
- Thick cushion base with different seating options depending on patient
- Tough, soft touch Dartex®
 Endurance material, with minimal stitching
- Strong, padded removable armrests

Specifications

| _ayflat dimensions: | | <u>1830mm x 660mm</u> | |
|---------------------|----------------------|-----------------------|--|
| | Seat width | up to 660mm | |
| | Height | 520-820mm | |
| | Hi-lo range | 490-790mm | |
| Wheel cas | /heel castors: total | | |
| & directional lock | | 125mm | |
| SWL | | 250 kg | |
| | | | |

Contact us for a demo or purchase:

updates@caremed.healthcare 0333 242 1250