


## ThermaZone® Pads

ThermaZone pads feature **Motion Advantage™** technology—leveraging precise fluid movement and delivering uniform temperature distribution across the surface of the pad.



FEATURING



Temperature controlled water  
flowing through welded  
microchannels

## Advantages of ThermaZone®

- Operates without the use of ice.
- Lightweight and portable
- 3 minutes to set up, requires only 2 oz. of distilled water.
- Operates at a constant temperature - has been tested for over 500 continuous hours
- Rapid achievement of temperature levels
- Pads that anatomically conform to all key muscle and joint areas
- Easy to setup, easy to use, easy to maintain
- Five preset timer options:
  - 10, 20 or 30 minutes
  - 30 minutes on and 30 minutes pause or
  - 30 minutes on and 60 minutes pause



**ThermaZone®**  
**Thermal Therapy**  
An opioid-free alternative




### Innovative Medical Equipment, LLC

23950 Commerce Park, Suite 100 | Cleveland, OH | 44122 USA  
(877) 901-ZONE (9663)



[www.therma-zone.com](http://www.therma-zone.com)

# M007 v8 04/2021



**heating + cooling therapy**  
when and where you need it most



**iceless. programmable. opioid-free.**



[www.therma-zone.com](http://www.therma-zone.com)

## Compelling Features


- Non-opioid therapy
- Heat or cool therapy in one device
- 3 minutes to set up and is virtually maintenance-free
- Short mode-changing delay enables contrast therapy at any time
- 20 adjustable comfort level settings

THERMAZONE COOL LEVEL										Temp. 72° F	THERMAZONE HEAT LEVEL									
10	9	8	7	6	5	4	3	2	1		1	2	3	4	5	6	7	8	9	10
34°	37°	40°	43°	46°	50°	54°	58°	62°	66°		80°	85°	90°	95°	100°	105°	110°	115°	120°	125°

Temperature delivery is based on an ambient temperature of 72° F.



**The Ice Age is Over!**  
No more ice, no more hydroculators.



**Broad Temp Range**  
Heats and cools; from 34°–125°F.  
20 comfort level settings.




**Smaller by Design**  
Lighter and smaller for easy portability.



**Motion Advantage Pads**  
For targeted therapy.



**Consistent Temp**  
Maintains temperature as long as you need it.



**H<sub>2</sub>O = 2 ounces**  
Uses only a 1/4 of a typical glass of water.



**Programmable Therapy**  
5 preset timer modes;  
heating and cooling therapy all in one device.



**Unit Size: 11" L x 4.5" W x 4.5" H**

**Unit Weight: 3 lbs. 7 oz.**



## Thermal Therapy Pads from Head to Toe

Individual pads work in conjunction with the ThermoZone® Thermal Therapy technology and are interchangeable with the device. All pads have adjustable straps for an excellent fit and light compression.

### Targeted Pads (Neoprene)



Knee, Elbow Pad



Back, Abdomen, Hip Pad



Range of Motion (ROM) Knee Brace



Shoulder Pad - Up to 44" Chest  
Large Shoulder Pad - Up to 54" Chest




Ankle Pad



Universal Pad

### General Pads



Front & Side Head Pad



Eye Pad



Occipital Pad  
(Back of Head)

### Accessories



Carrying Case  
(Holds Unit and 4 Pads)



Bed Bracket



IV Pole and Bed Bracket

### Single Patient Use Pads



Knee Pad



Back Pad



Large Shoulder Pad



Ankle Pad



Hip Pad