

for **orthopedics**

Featuring:  
**Single Patient Use Pads**

The **ThermaZone<sup>®</sup> Thermal Therapy Unit** and Single Patient Use Pads are an efficient combination to administer thermal therapy.

**ThermaZone** can be applied in the recovery room and is portable enough to take home for continued use.

**ThermaZone** does not require ice. It operates on two ounces of distilled water and is virtually maintenance-free.

**ThermaZone** has a broad temperature range from 34°F - 125°F, 20 comfort levels, and 5 pre-set timer modes.

Administering cold therapy post-surgery just got easier.



**ThermaZone<sup>®</sup>**

iceless. programmable. opioid-free.

**ThermaZone** is a powerful and compact device that uses thermoelectric technology to provide heating and cooling therapy where and when you need it, without the use of ice.

**ThermaZone** satisfies a clear unmet clinical need – targeted, localized, non-opioid therapy with consistent temperature management.



FEATURING




Temperature controlled water flowing through welded microchannels

Our pads feature **Motion Advantage<sup>™</sup> technology**, which leverages precise fluid movement to deliver uniform temperature distribution across the surface of the pad for maximum coverage.



Unit Size: 11" L x 4.5" W x 4.5" H  
Unit Weight: 3 lbs. 7 oz.




## Compelling Features

- Non-opioid therapy
- Heat or cool therapy in one device
- 3 minutes to set up requires only 2 oz. of distilled water
- Rapid achievement of temperature – 20 different comfort levels
- 5 preset timer options:
  - 10, 20, or 30 minutes
  - 30 minutes on and 30 minutes pause, or
  - 30 minutes on and 60 minutes pause



**The Ice Age Is Over!**  
No more ice, no more hydroculators.




**H<sub>2</sub>O = 2 ounces**  
Uses only 1/4 of a typical glass of water



**Motion Advantage Pads**  
Targeted therapy designed by surgeons.



**Smaller by Design**  
Lighter and smaller for easy portability.



**Broad Temp. Range**  
20 controlled temp settings from 34°–125°F.



**Programmable Therapy**  
5 preset timer modes; heating and cooling therapy all in one device



**Consistent Temp.**  
20 temp levels; maintains temp as long as you need it.

## Single Patient Use Pads (Disposable)

**ThermaZone Single Patient Use Pads** have the water delivery system welded directly into the pad for the perfect fit. They are durable, yet disposable. Pad interior is coated plastic which helps prevent fluid transfer.



**003-30**  
Shoulder Pad  
Fits up to a 54" Chest



**003-28**  
Back Pad



**003-27**  
Knee Pad



**003-29**  
Ankle Pad



Includes Unit, Power Supply, and Tubing Assembly

**ThermaZone Therapy System**

**003-99**  
Unit with Medical Grade Power Supply

## ThermaZone® Pads (Reusable)

Individual pads work in conjunction with the **ThermaZone Thermal Therapy** technology and are interchangeable with the unit. Pads are universal and fit right or left side.



**003-10**  
Front & Side Head Pad



**003-12**  
Eye Pad



**003-11**  
Occipital (Back of Head) Pad



**003-17**  
Knee | Elbow Pad



**003-47**  
Range of Motion (ROM) Knee Brace\*  
\*Insurance Reimbursable  
HPCS code L1833



**003-15**  
Shoulder Pad  
Fits up to a 44" Chest



**003-18**  
Back, Abdomen, Hip Pad




**003-22**  
Universal Pad



**003-19**  
Ankle | 90° Elbow Pad

**003-20**  
Large Shoulder Pad  
Fits up to a 54" Chest



**Innovative Medical Equipment, LLC**  
23950 Commerce Park, Suite 100 | Beachwood, OH 44122 USA  
(877) 901-ZONE (9663) | [www.therma-zone.com](http://www.therma-zone.com)

© 2021 Innovative Medical Equipment, LLC



#M003 v6 05/2021