ThermaZone Range of Motion (ROM) Knee Brace with Built-In Cooling



Introducing the **ThermaZone ROM Knee Brace** with an embedded water pathway that provides controlled temperature of both heating and cooling therapy. When used in conjunction with the ThermaZone Thermal Therapy device (ILO, 890.5720), patients benefit from both thermal therapy at controlled temperature and limited range of motion.

The ROM Knee Brace provides stability and support utilizing two polycentric hinges.

SPECIFICATIONS:

- Removable pins for angle adjustment:
 - ° Flexion Angles 0°, 10°, 20°, 40°, 60°, 90°
 - ° Extension Angles 0°, 10°, 20°, 40°, 60°, 90°, 150°
- Adjustable for various patient sizes (up to 26", 3" above the patella)

ThermaZone is a powerful and compact device that uses thermoelectric technology to provide heating and cooling therapy where and when you need it, without the use of ice.

ThermaZone satisfies a clear unmet clinical need - targeted, localized, non-opiod therapy with consistent temperature management.



Our pads feature Motion Advantage[™] technology, which leverages precise fluid movement to deliver uniform temperature distributions across the surface of the pad for maximum coverage.



PART #003-47 FEATURES:

- Insurance Reimbursable under HCPCS code L1833
- Made with durable neoprene material to provide support & compression
- Open patella design
- Fits left or right knee
- Comfortable & Lightweight
- Knee Brace with two polycentric hinges with adjustable flexion and extension joints
- A third strap behind the knee helps to keep the brace in place, requiring little to no adjustments during use
- Designed for patients who can bear weight on the knee



Unit Size: 11" L x 4.5" W x 4.5" H Unit Weight: 3lbs. 7 oz.





Compelling Features

- Non-opioid therapy
- · Heat or cool therapy in one device
- 3 minutes to set up, requires only 2 oz. of distilled water
- Rapid achievement of temperature 20 different comfort levels
- 5 preset timer options:
 - 10, 20, 30 minutes
- 30 minutes on and 30 minutes pause, or
- 30 minutes on and 60 minutes pause



The Ice Age is Over! No more ice.





Consistent Temp Maintains temperature as long as you need it.



Motion Advantage Pads Targeted therapy Designed by Surgeons.

ThermaZone Single Patient Use Pads have the water delivery system welded directly into the pad for the perfect fit. They are durable, yet disposable. Pad

Smaller by Design Lighter and smaller for easy portability.

Single Patient Use Pads





Pads Made with

ThermaZone®

003-15 Shoulder Pad - Fits up to a 44" Chest*

003-20 Large Shoulder Pad - Fits up to a 54" Chest*

003-28 Back Pad*

interior is coated with plastic which helps prevent fluid transfer. 003-29 Ankle Pad*

Individual pads work in conjunction with the ThermaZone Thermal Therapy technology and are interchangeable with the unit.





003-30 Large Shoulder Pad*

ThermaZone® Unit

Includes Unit and Tubing Assembly

ThermaZone Therapy System

003-99 - Unit with Medical Grade Power Supply





003-18 Back, Abdomen, Hip Pad*







*



003-11 Occipital (Back of Head) Pad





003-22 Universal Pad*



003-47 Range of Motion Knee Brace* with Polycentric Hinges

M025 v1 12/2020

Innovative Medical Equipment, LLC

23950 Commerce Park, Suite 100 • Cleveland, OH 44122 USA (877) 901-ZONE (9663)

003-19 Ankle Pad*















www.therma-zone.com ThermaZone is a registered trademark of

Innovative Medical Equipment, LLC © 2012

