

TrueStim

TREAT · RECOVER · PREVENT



What is TrueStim NMES?

TrueStim is an electric muscle stimulator utilized for treatment, recovery, and prevention of WC injuries. It is a device used for soft tissue healing, adjunctive relief of pain, swelling, and inflammation. It will insure healing more quickly and effectively, decreasing medication use, and allowing patient to return to full status more quickly. It has six channels with twelve treatment modes and is blue tooth wireless.

Treatment Indications

Muscle Weakness · Sore Joints · Muscle Spasticity · Prevention of Atrophy · Decreased Range of Motion · Pain Due to Muscle Weakness · Neuromuscular Deficits Effects Muscle Strength and Function

FEATURES

1

Combines microcurrent, Stim, and NMES into one machine.

2

Realistic traditional Chinese massage (Tuina) sensations and acupuncture technology.

3

Portable, lightweight, durable design.

4

Reusable, self-adhesive, silicon electrode gel pads.

5

Adjustable Timer.

DIRECTIONS



Place both electrode pads on the area of your body where you are experiencing stress, stiffness, aches, or pain.



Set the unit to the mode and intensity level that fits your personal preference.



Sit back and enjoy the deep, relaxing sensations of your TrueStim electronic muscle stimulator.

Benefits/Outcomes - Medical Necessity

Manage Inflammation · Reduce Muscle Spasms · Increase ROM · Reduce Joint Stiffness
Increase Circulation · Reduce Instability · Improve Muscle Functionality · Strengthen Muscles and Joints

TrueStim

FOR ORDERING INFORMATION, SPEAK WITH YOUR EASTERN ORTHOPAEDICS REPRESENTATIVE AT (800) 366-8051. ALSO VISIT OUR WEBSITE AT EORTHO.COM