


# **BIOSOLID KIT EXAMPLE**



**BioSolid 503 Complete + Mercury**



**Pathogens**



**PFAS**

## **Included in the BioSolids Kit are 10 bottles:**

503 Complete + Mercury: 1000 mL and 250 mL unpreserved bottles, both bottles must be used.

Pathogens: 7 sterile bottles with Sodium Thiosulfate.

PFAS: Sterile 250 mL bottle with Barcode label.

**ALL BOTTLES MUST BE FILLED  
FOR REQUESTED TESTING**

Please make sure all labels and Chain of Custodies are correctly filled out with Client information (Client, Date, Time, and Analytes).