

INNOVATIVE PREVENTATIVE CARE

The Capiflo NS9000 serves as a portable solution for Deep Vein Thrombosis (DVT) Prevention Therapy. It emulates the body's intrinsic muscle contractions, effectively amplifying blood circulation and augmenting venous blood flow to the extremities. This portable device ensures uninterrupted DVT prevention for more than 24 hours, facilitating patient mobility and independence. While it comes pre-set at 55 mmHg pressure, the NS9000 allows for personalized adjustments, catering to optimal patient well-being.

TECHNICAL DATA

-ŒD-Power Supply

Class II, input: 100-240 Vac 50-60 Hz Output 7.4 v dc @1.1 amps

Battery Charge

Lasts approximately 24 hours. Recharge requires 4 hours

Default Setting

55 mmHg
Alternate Settings: 50 mmHg
or 40 mmHg by depressing button

Controller unit 2.6" x 5.2"

Weight .5 lb sleeve

Mode of operation: Cyclical

Source of Power 7.4 Volt Li-ion Battery

2-3.7 volt cells