

# **Sili-Fert A**

*Silicon for race pigeons*

Choline stabilized  
Orthosilicic Acid



For healthy chicks  
and strong pigeons

## **PIGEONS**

TECHNOLOGICAL ADDITIVE  
animal food supplement

**Sili-Fert A**  
*Silicon for race pigeons*



Choline stabilized  
Orthosilicic Acid



For the hobby pigeon fancier  
and for professional pigeon sport

**Professional**  
Silicon Food Supplement



GLUCOSAMINE

COLLAGEN



Ch-OSA Choline Orthosilicic Acid  
[www.liquidsilicium.com](http://www.liquidsilicium.com)  
promotes the biosynthesis of...

CREATINE

**Sili-Fert A**  
Silicon for race pigeons



IMMUNITY

DIGESTION



LESS  
ANTIBIOTICS

**Sili-Fert A**  
— Silicon for race pigeons



**Sili-Fert A**  
Silicon for race pigeons



Ch-OSA Choline stabilized  
Ortho-Silicic Acid  
2.8% Si [OH]<sub>4</sub>

Liquid Silicium : activates + regenerates + stimulates + regulates + protects

- ✓ maintenance of the joints, tendons, muscles, bones...
- ✓ improves the absorption of minerals and vitamins
- ✓ promotes the digestions, metabolism and growth
- ✓ prevents weak and translucent animals
- ✓ speeds up the healing proces

**Choline stabilized Orthosilicic Acid**

[www.liquidsilicium.com](http://www.liquidsilicium.com)

Professional nutrient for  
**PIGEONS**

**Sili-Fert A**  
*Silicon for race pigeons*

STRONG  
BONES

[www.liquidilicium.com](http://www.liquidilicium.com)



OPTIMAL  
RESISTANCE

FAST  
RECOVERY

# Sili-Fert A

Silicon for race pigeons

Choline stabilized  
Orthosilicic Acid

## Daily dosis

1 to 2 drops per 5 kg

**34,99 €**

(e: 100 ml)

## Dilution

1 ml per 1000ml (1/1000)

or 20 - 30 drops per liter

opm: 20 drops = 1 ml

## Administration

Dilute Sili-Fert A in the drinking water or dilute in a little water and sprinkle Sili-Fert A daily on the feed



Admission: BE355

EAN: 6090950523524

**PIGEON**  
Technological Additive



info@silifert.eu  
+32 (0)468 115 335



# Sili-Fert A

Silicon for race pigeons

Choline stabilized  
Orthosilicic Acid

Liquid Silicium : activates + regenerates + stimulates + regulates + protects



## Dilution

Add daily one drop per 5 kg or 10 pounds bodymass

- dilute Sili-Fert A in a little water and mix it in the feed  
or
- dilute Sili-Fert A into the drinking water @ 1/1000 :  
e.g. 100 liter water = 100.000 ml -> 1/1000 = 100 ml Sili-Fert A

## Administration

(20 drops = 1 ml)

Estimate the total weight of the chickens (#1 drop per 5 kg),  
dilute Sili-Fert P in water and sprinkle it on the animal feed

1. ex. 20 kg or 40 pounds = 4 drops, diluted in minimum 10 ml
2. ex. 100 kg (200 pounds) = 20 drops (= 1 ml) diluted in 1000 ml  
(20 drops = 1 ml)



# Sili-Fert A

Silicon for race pigeons

Choline stabilized

OrthoSilicic Acid



100 ml dropper - 1000 ml bottle \* 5 - 10 - 20 liter can \* 60 liter drum \* 600 - 800 - 1000 L IBC



ing. Robin Severi

robin@severi.technology

+32 468 115 335

International distributor  
Severi Technologies - AquaNano  
info@silifert.eu

[www.liquidsilicium.com](http://www.liquidsilicium.com)

Liquid Silicium : 100% Bio-Available Silicon Technological Nutrition Additive for Pet Animals