

Belmont^{tb}Compass with Delivery

Left / Right Convertible Treatment Centre

100% Japanese Made
**WITH FIVE YEAR
EQUIPMENT WARRANTY**

Belmont

Reliability - Think Ahead

The tbCompass' flexible systems offer surgeries the versatility of an ambidextrous Operator's console or cuspidor. Its Left-Right conversion pivot at the head of the chair requires less space in the surgery when switching from Left to Right. It also has flexible configurations to suit any practice and procedure thus allowing it to fit in smaller spaces. On top of that, it offers flexible instrument configurations to fulfill individual preferences.

- Ambidextrous: Left or Right handed set up
- Air and Electric Operating System
- Seamless or Plush Upholstery option
- 10 LED Dental Light with touchless on/off sensor
- Twin-axis Headrest and Ultra Slim Backrest to improve operator's working comfort
- Ambidextrous Cuspidor with rotating and detachable spittoon bowl

60 Second Changeover

From left to right-side, or vice versa requires no tools or expertise. The simple procedure results in identical conditions for left and right handed operators.




Right-handed set up.
Raise chair height and remove the rear cover.



Incline chair.
1-2. Swing table round to left side.



3-4. Swing cuspidor round to right side.



5. Secure cuspidor in position.
Secure cuspidor locking lever.
2. Rotate console.
Left-handed set up.



* Handpieces shown are for illustration purposes only, we will customise specifically to your needs. Ask us about optional extras and package deals.