


















ALPRO

Product Image	Product Code	Product Detail	NZD RRP + GST
CLEANING & DISINFECTION OF ASPIRATION SYSTEMS			
	CH ALPRO-3101	<p align="center">AlproJet-D 5 Litre</p> <p>Highly efficient, high pH, aldehyde-free and non-foaming liquid concentrate for the daily cleaning and maintenance of dental aspiration and separation systems, with or without amalgam separator</p>	\$129.95
	CH ALPRO-3105	<p align="center">AlproJet-W 5 Litre</p> <p>Highly efficient, low pH, aldehyde-free and non-foaming liquid concentrate for the weekly cleaning and maintenance of dental aspiration and separation systems, with and without amalgam separator.</p>	\$139.95
	CH ALPRO-3104	<p align="center">AlproJet-D 1 Litre</p> <p>Highly efficient, high pH, aldehyde-free and non-foaming liquid concentrate for the daily cleaning and maintenance of dental aspiration and separation systems, with or without amalgam separator</p>	\$51.95
	CH ALPRO-3106	<p align="center">AlproJet-W 1 Litre</p> <p>Highly efficient, low pH, aldehyde-free and non-foaming liquid concentrate for the weekly cleaning and maintenance of dental aspiration and separation systems, with and without amalgam separator.</p>	\$40.95

	<p>CH ALPRO-4658</p>	<p>IC-100 500ml Spray Bottle Neutral detergent for all surfaces and materials. Non protein fixing neutral detergent and cleaning foam, ready for use for the cleaning of hard surfaces, dental units, x-ray machines, all turbines as well as pre-cleaning of intraoral instruments. Dissolves residues of blood and other body fluids from surfaces of medical devices and instruments. Stops protein fixation on instruments before ultrasonic cleaning.</p>	<p>TBC</p>
	<p>CH ALPRO-4659</p>	<p>IC-100 5 Litre Canister Neutral detergent for all surfaces and materials. Non protein fixing neutral detergent and cleaning foam, ready for use for the cleaning of hard surfaces, dental units, x-ray machines, all turbines as well as pre-cleaning of intraoral instruments. Dissolves residues of blood and other body fluids from surfaces of medical devices and instruments. Stops protein fixation on instruments before ultrasonic cleaning</p>	<p>\$123.95</p>
	<p>CH ALPRO-4625</p>	<p>PlastiSept eco tub of 70 wipes Ready to use alcohol-free cleaning and disinfecting solution without any aldehydes, phenols and QAC. Especially suitable for alcohol sensitive surfaces of medical devices such as dental units, suspension tables, x-ray machines, counter tops, hand pieces and turbines.</p>	<p>\$49.95</p>
	<p>CH ALPRO-4624</p>	<p>PlastiSept eco refills 70 wipes</p>	<p>\$38.95</p>

	<p>CH ALPRO-3024</p>	<p align="center">BIB forte – 1 Litre</p> <p align="center">Instrument Grade Disinfectant – Liquid concentrate free of any aldehydes and phenols, designed for non-protein fixing cleaning and disinfection of non-critical instruments and medical devices.</p>	<p align="center">\$129.95</p>
	<p>CH ALPRO-3183</p>	<p align="center">Alpron 1 Litre</p> <p align="center">Liquid concentrate for decontamination of procedural water in dental treatment units.</p>	<p align="center">\$129.95</p>
	<p>CH ALPRO-3181</p>	<p align="center">Bilpron 1 Litre</p> <p align="center">Biodegradable, ready to use solution for the decontamination and prevention of formation of biofilms in procedural water lines in medical and dental treatment units.</p>	<p align="center">\$55.95</p>
	<p>CH ALPRO-3179</p>	<p align="center">BRS Bilpron Starter Set</p> <p align="center">Biofilm Removal System includes Bilpron, BRS PreCleaner and BRS Remover</p>	<p align="center">\$175.95</p>

	<p>CH ALPRO-3199</p>	<p align="center">BC-San 100ml 12 Pack 12 x BC-San 100ml bottles, cleaning of the water bottle and GBT cylinder</p>	<p align="right">\$16.95</p>
	<p>CH ALPRO-4153</p>	<p align="center">CleanWipes – 70 Wipes – Cleaning of medical devices Alcohol-free wet wipes for the non-protein fixing cleaning and decontamination of surfaces of medical devices. CleanWipes remove easily and gently organic and inorganic soilings. CleanWipes (soaked with IC-100) solution are ready for use cleaning wipes, especially designed for cleaning and decontamination.</p>	<p align="right">\$49.95</p>
	<p>CH ALPRO-4154</p>	<p align="center">CleanWipes IC-100 refills 70 wipes</p>	<p align="right">\$28.95</p>
	<p>CH ALPRO-3079</p>	<p align="center">Biotesters ea.</p>	<p align="right">\$29.95</p>
	<p>CH ALPRO-3085</p>	<p align="center">Multifunctional Container Multifunctional Container for the solutions used for cleaning and maintenance of the aspiration system.</p>	<p align="right">\$59.95</p>

