



Adhesive Code	Filler Type	Number of Components	Work Life @ 25°C	Cure Schedule	Volume Resistivity	Extractable Chloride Ionic Content	Die Shear Strength	Price
ZTERIAL 81	Clear	2	2 Hours	120°C, 5 Minutes		< 15 ppm	> 10 Kg	\$0.32/g
ZTERIAL 82	Clear	2	3 Hours	120°C, 5 Minutes		< 15 ppm	> 10 Kg	\$0.32/g
ZTERIAL 91	Non-Silver	2	3 Hours	120°C, 10 Minutes		< 15 ppm	> 10 Kg	\$0.43/g
ZTERIAL 92	Non-Silver	1	14 Days	150°C, 1 Hours		< 15 ppm	> 10 Kg	\$0.43/g
ZTERIAL 122	Silver	2	3 Hours	120°C, 5 Minutes	< 0.00007 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 317	Silver	2	8 Hours	120°C, 5 Minutes	< 0.00006 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 322	Silver	2	8 Hours	120°C, 5 Minutes	< 0.00007 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 420	Silver	2	18 Hours	120°C, 15 Minutes	< 0.0003 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 421	Silver	2	18 Hours	120°C, 15 Minutes	< 0.0001 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 422	Silver	2	18 Hours	120°C, 15 Minutes	< 0.00007 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 515	Silver	2	3 Days	150°C, 1 Hours	< 0.00008 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 516	Silver	2	3 Days	150°C, 1 Hours	< 0.00006 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 517	Silver	2	3 Days	150°C, 1 Hours	< 0.00006 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 518	Silver	2	3 Days	150°C, 1 Hours	< 0.0003 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 519	Silver	2	3 Days	150°C, 1 Hours	< 0.0005 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 615	Silver	1	14 Days	150°C, 1 Hours	< 0.00006 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 616	Silver	1	14 Days	150°C, 1 Hours	< 0.00008 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 617	Silver	1	14 Days	150°C, 1 Hours	< 0.0001 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 618	Silver	1	14 Days	150°C, 1 Hours	< 0.00008 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g
ZTERIAL 619	Silver	1	14 Days	150°C, 1 Hours	< 0.0005 ohms-cm	< 15 ppm	> 10 Kg	\$1.62/g